

Table of contents

\\QIS_LAB

BME Lab

2026

Group 4

MEGRE_FA9_TR12_TE1_1p33_TE2_3p27
MEGRE_FA9_TR14_TE1_2p41_TE2_5p42
MEGRE_FA20_TR14_TE1_2p41_TE2_5p42
se_multiecho_TR1500ms
se_multiecho_TR5000ms

\\QIS_LAB\BME Lab\2026\Group 4\MEGRE_FA9_TR12_TE1_1p33_TE2_3p27

TA: 1:35 PM: FIX Voxel size: 1.4×1.4×3.0 mmPAT: Off Rel. SNR: 1.00 : tfi

Properties

Prio recon	Off
Load images to viewer	On
Inline movie	Off
Auto store images	On
Load images to stamp segments	Off
Load images to graphic segments	Off
Auto open inline display	Off
Auto close inline display	Off
Start measurement without further preparation	Off
Wait for user to start	Off
Start measurements	Single measurement

Routine

Slab group	1
Slabs	1
Dist. factor	20 %
Position	L1.7 P59.0 H23.3 mm
Orientation	Coronal
Phase enc. dir.	R >> L
AutoAlign	Head > Brain
Phase oversampling	0 %
Slice oversampling	0.0 %
Slices per slab	48
FoV read	224 mm
FoV phase	100.0 %
Slice thickness	3.00000 mm
TR	1928.00 ms
TE 1	1.33 ms
TE 2	3.27 ms
TE 3	5.14 ms
TE 4	7.01 ms
TE 5	8.88 ms
TE 6	10.75 ms
Averages	1
Concatenations	1
Filter	None
Coil elements	BO1-3;SP2-4

Contrast - Common

TR	1928.00 ms
TE 1	1.33 ms
TE 2	3.27 ms
TE 3	5.14 ms
TE 4	7.01 ms
TE 5	8.88 ms
TE 6	10.75 ms
Magn. preparation	None
Flip angle	9 deg
Fat suppr.	None
Wrap-up Magn.	None

Contrast - Dynamic

Averages	1
Averaging mode	Short term
Reconstruction	Magn./Phase
Measurements	1
Multiple series	Off

Resolution - Common

FoV read	224 mm
FoV phase	100.0 %
Slice thickness	3.00000 mm
Base resolution	160
Phase resolution	100 %
Slice resolution	100 %
Phase partial Fourier	Off
Slice partial Fourier	Off
Trajectory	Cartesian
Interpolation	Off

Resolution - iPAT

PAT mode	None
----------	------

Resolution - Filter Image

Image Filter	Off
Distortion Corr.	Off
Prescan Normalize	Off
Normalize	Off
B1 filter	Off

Resolution - Filter Rawdata

Raw filter	Off
Elliptical filter	Off
POCS	Off

Geometry - Common

Slab group	1
Slabs	1
Dist. factor	20 %
Position	L1.7 P59.0 H23.3 mm
Orientation	Coronal
Phase enc. dir.	R >> L
Slice oversampling	0.0 %
Slices per slab	48
FoV read	224 mm
FoV phase	100.0 %
Slice thickness	3.00000 mm
TR	1928.00 ms
Multi-slice mode	Sequential
Series	Interleaved
Concatenations	1

Geometry - AutoAlign

Slab group	1
Position	L1.7 P59.0 H23.3 mm
Orientation	Coronal
Phase enc. dir.	R >> L
AutoAlign	Head > Brain
Initial Position	L1.7 P59.0 F12.7
L	1.7 mm
P	59.0 mm
H	12.7 mm
Initial Rotation	0.00 deg
Initial Orientation	Coronal

Geometry - Saturation

Fat suppr.	None
Wrap-up Magn.	None
Special sat.	None

Geometry - Navigator**Geometry - Tim Planning Suite**

Set-n-Go Protocol	Off
Table position	H
Table position	36 mm
Inline Composing	Off

System - Miscellaneous

Positioning mode	FIX
Table position	H
Table position	36 mm
MSMA	S - C - T
Sagittal	R >> L
Coronal	A >> P
Transversal	F >> H
Coil Combine Mode	Adaptive Combine
Save uncombined	Off
Matrix Optimization	Off
Coil Focus	Flat
AutoAlign	Head > Brain
Coil Select Mode	Default

System - Adjustments

B0 Shim mode	Standard
B1 Shim mode	Patient-specific
Adjust with body coil	Off
Confirm freq. adjustment	Off
Assume Dominant Fat	Off
Assume Silicone	Off
Adjustment Tolerance	Auto

System - Adjust Volume

! Position	L1.5 P72.5 H13.0 mm
! Orientation	S > T1.0 > C0.2
! Rotation	-0.94 deg
! A >> P	44 mm
! F >> H	137 mm
! R >> L	87 mm
Reset	Off

System - pTx Volumes

B1 Shim mode	Patient-specific
Excitation	Non-sel.

System - Tx/Rx

Frequency 1H	123.229536 MHz
Correction factor	1
Gain	Low
Img. Scale Cor.	1.000
Reset	Off
? Ref. amplitude 1H	0.000 V

Physio - Signal1

1st Signal/Mode	None
TR	1928.00 ms
Concatenations	1
Segments	160

Physio - Cardiac

Tagging	None
Magn. preparation	None
Fat suppr.	None
Dark blood	Off

Physio - Cardiac

FoV read	224 mm
FoV phase	100.0 %
Phase resolution	100 %
Cine	Off
Trajectory	Cartesian
Dummy heartbeats	1

Physio - PACE

Resp. control	Off
Concatenations	1

Inline - Common

Subtract	Off
Measurements	1
StdDev	Off
Save original images	On

Inline - Cardiac

Inline Evaluation	Off
Magn. preparation	None
Contrasts	6
TE 1	1.33 ms
TE 2	3.27 ms
TE 3	5.14 ms
TE 4	7.01 ms
TE 5	8.88 ms
TE 6	10.75 ms
TR	1928.00 ms
Save original images	On

Inline - MIP

MIP-Sag	Off
MIP-Cor	Off
MIP-Tra	Off
MIP-Time	Off
Save original images	On

Inline - Composing

Inline Composing	Off
Distortion Corr.	Off

Sequence - Part 1

Introduction	Off
Dimension	3D
Elliptical scanning	Off
Reordering	Linear
Asymmetric echo	Off
Contrasts	6
Flow comp. 1	No
Readout mode	Monopolar
Optimization	Min. TE TR
Multi-slice mode	Sequential
Echo spacing	12.1 ms
Sequence type	Gre
Bandwidth 1	1008 Hz/Px
Bandwidth 2	1953 Hz/Px
Bandwidth 3	1953 Hz/Px
Bandwidth 4	1953 Hz/Px
Bandwidth 5	1953 Hz/Px
Bandwidth 6	1953 Hz/Px

Sequence - Part 2

Define	Segments
--------	----------

Sequence - Part 2

Segments	160
RF pulse type	Fast
Gradient mode	Performance
Excitation	Non-sel.
Flip angle mode	Constant
RF spoiling	On
Phase Enc. Rewinder	On
Cine	Off

Sequence - Special

3D Rad Acq. [L] # 0	1 n.a.
3D Rad Acq. [L] # 1	2 n.a.
3D Rad Acq. [L] # 2	2 n.a.
3D Rad Rec. [L] # 0	2 n.a.
3D Rad Rec. [L] # 1	4 n.a.
3D Rad Rec. [L] # 2	2 n.a.
3D Rad Rec. [L] # 3	2 n.a.
3D Rad Rec. [L] # 4	75 n.a.
3D Rad Rec. [D] # 0	1.0 n.a.
3D Rad Rec. [D] # 1	10.0 n.a.
Prep Pulses [L] # 0	0 n.a.
LIBRE Pulse [L] # 0	1 n.a.
LIBRE Pulse [L] # 1	0 n.a.
LIBRE Pulse [L] # 2	300 n.a.
LIBRE Pulse [L] # 3	200 n.a.
LIBRE Pulse [L] # 4	0 n.a.
LIBRE Pulse [L] # 5	0 n.a.
LIBRE Pulse [L] # 6	260 n.a.
Self-Navigation [L] # 0	2 n.a.
Self-Navigation [L] # 1	2 n.a.
Self-Navigation [L] # 2	1 n.a.
Self-Navigation [L] # 3	0 n.a.
Self-Navigation [L] # 4	2 n.a.
Self-Navigation [L] # 5	2 n.a.
Self-Navigation [L] # 6	80 n.a.
Self-Navigation [L] # 7	20 n.a.
Self-Navigation [L] # 8	20 n.a.
Self-Navigation [L] # 9	1 n.a.
Self-Gating [L] 0	2 n.a.

Sequence - Assistant

Mode	Off
Allowed delay	0 s

\\QIS_LAB\BME Lab\2026\Group 4\MEGRE_FA9_TR14_TE1_2p41_TE2_5p42

TA: 1:53 PM: FIX Voxel size: 1.4×1.4×3.0 mmPAT: Off Rel. SNR: 1.00 : tfi

Properties

Prio recon	Off
Load images to viewer	On
Inline movie	Off
Auto store images	On
Load images to stamp segments	Off
Load images to graphic segments	Off
Auto open inline display	Off
Auto close inline display	Off
Start measurement without further preparation	Off
Wait for user to start	Off
Start measurements	Single measurement

Routine

Slab group	1
Slabs	1
Dist. factor	20 %
Position	L1.7 P59.0 H23.3 mm
Orientation	Coronal
Phase enc. dir.	R >> L
AutoAlign	Head > Brain
Phase oversampling	0 %
Slice oversampling	0.0 %
Slices per slab	48
FoV read	224 mm
FoV phase	100.0 %
Slice thickness	3.00000 mm
TR	2299.20 ms
TE 1	2.41 ms
TE 2	5.42 ms
TE 3	7.33 ms
TE 4	9.24 ms
TE 5	11.15 ms
TE 6	13.06 ms
Averages	1
Concatenations	1
Filter	None
Coil elements	BO1-3;SP2-4

Contrast - Common

TR	2299.20 ms
TE 1	2.41 ms
TE 2	5.42 ms
TE 3	7.33 ms
TE 4	9.24 ms
TE 5	11.15 ms
TE 6	13.06 ms
Magn. preparation	None
Flip angle	9 deg
Fat suppr.	None
Wrap-up Magn.	None

Contrast - Dynamic

Averages	1
Averaging mode	Short term
Reconstruction	Magn./Phase
Measurements	1
Multiple series	Off

Resolution - Common

FoV read	224 mm
FoV phase	100.0 %
Slice thickness	3.00000 mm
Base resolution	160
Phase resolution	100 %
Slice resolution	100 %
Phase partial Fourier	Off
Slice partial Fourier	Off
Trajectory	Cartesian
Interpolation	Off

Resolution - iPAT

PAT mode	None
----------	------

Resolution - Filter Image

Image Filter	Off
Distortion Corr.	Off
Prescan Normalize	Off
Normalize	Off
B1 filter	Off

Resolution - Filter Rawdata

Raw filter	Off
Elliptical filter	Off
POCS	Off

Geometry - Common

Slab group	1
Slabs	1
Dist. factor	20 %
Position	L1.7 P59.0 H23.3 mm
Orientation	Coronal
Phase enc. dir.	R >> L
Slice oversampling	0.0 %
Slices per slab	48
FoV read	224 mm
FoV phase	100.0 %
Slice thickness	3.00000 mm
TR	2299.20 ms
Multi-slice mode	Sequential
Series	Interleaved
Concatenations	1

Geometry - AutoAlign

Slab group	1
Position	L1.7 P59.0 H23.3 mm
Orientation	Coronal
Phase enc. dir.	R >> L
AutoAlign	Head > Brain
Initial Position	L1.7 P59.0 F12.7
L	1.7 mm
P	59.0 mm
H	12.7 mm
Initial Rotation	0.00 deg
Initial Orientation	Coronal

Geometry - Saturation

Fat suppr.	None
Wrap-up Magn.	None
Special sat.	None

Geometry - Navigator**Geometry - Tim Planning Suite**

Set-n-Go Protocol	Off
Table position	H
Table position	36 mm
Inline Composing	Off

System - Miscellaneous

Positioning mode	FIX
Table position	H
Table position	36 mm
MSMA	S - C - T
Sagittal	R >> L
Coronal	A >> P
Transversal	F >> H
Coil Combine Mode	Adaptive Combine
Save uncombined	Off
Matrix Optimization	Off
Coil Focus	Flat
AutoAlign	Head > Brain
Coil Select Mode	Default

System - Adjustments

B0 Shim mode	Standard
B1 Shim mode	Patient-specific
Adjust with body coil	Off
Confirm freq. adjustment	Off
Assume Dominant Fat	Off
Assume Silicone	Off
Adjustment Tolerance	Auto

System - Adjust Volume

! Position	L1.5 P72.5 H13.0 mm
! Orientation	S > T1.0 > C0.2
! Rotation	-0.94 deg
! A >> P	44 mm
! F >> H	137 mm
! R >> L	87 mm
Reset	Off

System - pTx Volumes

B1 Shim mode	Patient-specific
Excitation	Non-sel.

System - Tx/Rx

Frequency 1H	123.229536 MHz
Correction factor	1
Gain	Low
Img. Scale Cor.	1.000
Reset	Off
? Ref. amplitude 1H	0.000 V

Physio - Signal1

1st Signal/Mode	None
TR	2299.20 ms
Concatenations	1
Segments	160

Physio - Cardiac

Tagging	None
Magn. preparation	None
Fat suppr.	None
Dark blood	Off

Physio - Cardiac

FoV read	224 mm
FoV phase	100.0 %
Phase resolution	100 %
Cine	Off
Trajectory	Cartesian
Dummy heartbeats	1

Physio - PACE

Resp. control	Off
Concatenations	1

Inline - Common

Subtract	Off
Measurements	1
StdDev	Off
Save original images	On

Inline - Cardiac

Inline Evaluation	Off
Magn. preparation	None
Contrasts	6
TE 1	2.41 ms
TE 2	5.42 ms
TE 3	7.33 ms
TE 4	9.24 ms
TE 5	11.15 ms
TE 6	13.06 ms
TR	2299.20 ms
Save original images	On

Inline - MIP

MIP-Sag	Off
MIP-Cor	Off
MIP-Tra	Off
MIP-Time	Off
Save original images	On

Inline - Composing

Inline Composing	Off
Distortion Corr.	Off

Sequence - Part 1

Introduction	Off
Dimension	3D
Elliptical scanning	Off
Reordering	Linear
Asymmetric echo	Off
Contrasts	6
Flow comp. 1	No
Readout mode	Monopolar
Optimization	Min. TE TR
Multi-slice mode	Sequential
Echo spacing	14.4 ms
Sequence type	Gre
Bandwidth 1	300 Hz/Px
Bandwidth 2	1953 Hz/Px
Bandwidth 3	1953 Hz/Px
Bandwidth 4	1953 Hz/Px
Bandwidth 5	1953 Hz/Px
Bandwidth 6	1953 Hz/Px

Sequence - Part 2

Define	Segments
--------	----------

Sequence - Part 2

Segments	160
RF pulse type	Fast
Gradient mode	Performance*
Excitation	Non-sel.
Flip angle mode	Constant
RF spoiling	On
Phase Enc. Rewinder	On
Cine	Off

Sequence - Special

3D Rad Acq. [L] # 0	1 n.a.
3D Rad Acq. [L] # 1	2 n.a.
3D Rad Acq. [L] # 2	2 n.a.
3D Rad Rec. [L] # 0	2 n.a.
3D Rad Rec. [L] # 1	4 n.a.
3D Rad Rec. [L] # 2	2 n.a.
3D Rad Rec. [L] # 3	2 n.a.
3D Rad Rec. [L] # 4	75 n.a.
3D Rad Rec. [D] # 0	1.0 n.a.
3D Rad Rec. [D] # 1	10.0 n.a.
Prep Pulses [L] # 0	0 n.a.
LIBRE Pulse [L] # 0	1 n.a.
LIBRE Pulse [L] # 1	0 n.a.
LIBRE Pulse [L] # 2	300 n.a.
LIBRE Pulse [L] # 3	200 n.a.
LIBRE Pulse [L] # 4	0 n.a.
LIBRE Pulse [L] # 5	0 n.a.
LIBRE Pulse [L] # 6	260 n.a.
Self-Navigation [L] # 0	2 n.a.
Self-Navigation [L] # 1	2 n.a.
Self-Navigation [L] # 2	1 n.a.
Self-Navigation [L] # 3	0 n.a.
Self-Navigation [L] # 4	2 n.a.
Self-Navigation [L] # 5	2 n.a.
Self-Navigation [L] # 6	80 n.a.
Self-Navigation [L] # 7	20 n.a.
Self-Navigation [L] # 8	20 n.a.
Self-Navigation [L] # 9	1 n.a.
Self-Gating [L] 0	2 n.a.

Sequence - Assistant

Mode	Off
Allowed delay	0 s

\\QIS_LAB\BME Lab\2026\Group 4\MEGRE_FA20_TR14_TE1_2p41_TE2_5p42

TA: 1:53 PM: FIX Voxel size: 1.4×1.4×3.0 mmPAT: Off Rel. SNR: 1.00 : tfl

Properties

Prio recon	Off
Load images to viewer	On
Inline movie	Off
Auto store images	On
Load images to stamp segments	Off
Load images to graphic segments	Off
Auto open inline display	Off
Auto close inline display	Off
Start measurement without further preparation	Off
Wait for user to start	Off
Start measurements	Single measurement

Routine

Slab group	1
Slabs	1
Dist. factor	20 %
Position	L1.7 P59.0 H23.3 mm
Orientation	Coronal
Phase enc. dir.	R >> L
AutoAlign	Head > Brain
Phase oversampling	0 %
Slice oversampling	0.0 %
Slices per slab	48
FoV read	224 mm
FoV phase	100.0 %
Slice thickness	3.00000 mm
TR	2299.20 ms
TE 1	2.41 ms
TE 2	5.42 ms
TE 3	7.33 ms
TE 4	9.24 ms
TE 5	11.15 ms
TE 6	13.06 ms
Averages	1
Concatenations	1
Filter	None
Coil elements	BO1-3;SP2-4

Contrast - Common

TR	2299.20 ms
TE 1	2.41 ms
TE 2	5.42 ms
TE 3	7.33 ms
TE 4	9.24 ms
TE 5	11.15 ms
TE 6	13.06 ms
Magn. preparation	None
Flip angle	20 deg
Fat suppr.	None
Wrap-up Magn.	None

Contrast - Dynamic

Averages	1
Averaging mode	Short term
Reconstruction	Magn./Phase
Measurements	1
Multiple series	Off

Resolution - Common

FoV read	224 mm
FoV phase	100.0 %
Slice thickness	3.00000 mm
Base resolution	160
Phase resolution	100 %
Slice resolution	100 %
Phase partial Fourier	Off
Slice partial Fourier	Off
Trajectory	Cartesian
Interpolation	Off

Resolution - iPAT

PAT mode	None
----------	------

Resolution - Filter Image

Image Filter	Off
Distortion Corr.	Off
Prescan Normalize	Off
Normalize	Off
B1 filter	Off

Resolution - Filter Rawdata

Raw filter	Off
Elliptical filter	Off
POCS	Off

Geometry - Common

Slab group	1
Slabs	1
Dist. factor	20 %
Position	L1.7 P59.0 H23.3 mm
Orientation	Coronal
Phase enc. dir.	R >> L
Slice oversampling	0.0 %
Slices per slab	48
FoV read	224 mm
FoV phase	100.0 %
Slice thickness	3.00000 mm
TR	2299.20 ms
Multi-slice mode	Sequential
Series	Interleaved
Concatenations	1

Geometry - AutoAlign

Slab group	1
Position	L1.7 P59.0 H23.3 mm
Orientation	Coronal
Phase enc. dir.	R >> L
AutoAlign	Head > Brain
Initial Position	L1.7 P59.0 F12.7
L	1.7 mm
P	59.0 mm
H	12.7 mm
Initial Rotation	0.00 deg
Initial Orientation	Coronal

Geometry - Saturation

Fat suppr.	None
Wrap-up Magn.	None
Special sat.	None

Geometry - Navigator**Geometry - Tim Planning Suite**

Set-n-Go Protocol	Off
Table position	H
Table position	36 mm
Inline Composing	Off

System - Miscellaneous

Positioning mode	FIX
Table position	H
Table position	36 mm
MSMA	S - C - T
Sagittal	R >> L
Coronal	A >> P
Transversal	F >> H
Coil Combine Mode	Adaptive Combine
Save uncombined	Off
Matrix Optimization	Off
Coil Focus	Flat
AutoAlign	Head > Brain
Coil Select Mode	Default

System - Adjustments

B0 Shim mode	Standard
B1 Shim mode	Patient-specific
Adjust with body coil	Off
Confirm freq. adjustment	Off
Assume Dominant Fat	Off
Assume Silicone	Off
Adjustment Tolerance	Auto

System - Adjust Volume

! Position	L1.5 P72.5 H13.0 mm
! Orientation	S > T1.0 > C0.2
! Rotation	-0.94 deg
! A >> P	44 mm
! F >> H	137 mm
! R >> L	87 mm
Reset	Off

System - pTx Volumes

B1 Shim mode	Patient-specific
Excitation	Non-sel.

System - Tx/Rx

Frequency 1H	123.229536 MHz
Correction factor	1
Gain	Low
Img. Scale Cor.	1.000
Reset	Off
? Ref. amplitude 1H	0.000 V

Physio - Signal1

1st Signal/Mode	None
TR	2299.20 ms
Concatenations	1
Segments	160

Physio - Cardiac

Tagging	None
Magn. preparation	None
Fat suppr.	None
Dark blood	Off

Physio - Cardiac

FoV read	224 mm
FoV phase	100.0 %
Phase resolution	100 %
Cine	Off
Trajectory	Cartesian
Dummy heartbeats	1

Physio - PACE

Resp. control	Off
Concatenations	1

Inline - Common

Subtract	Off
Measurements	1
StdDev	Off
Save original images	On

Inline - Cardiac

Inline Evaluation	Off
Magn. preparation	None
Contrasts	6
TE 1	2.41 ms
TE 2	5.42 ms
TE 3	7.33 ms
TE 4	9.24 ms
TE 5	11.15 ms
TE 6	13.06 ms
TR	2299.20 ms
Save original images	On

Inline - MIP

MIP-Sag	Off
MIP-Cor	Off
MIP-Tra	Off
MIP-Time	Off
Save original images	On

Inline - Composing

Inline Composing	Off
Distortion Corr.	Off

Sequence - Part 1

Introduction	Off
Dimension	3D
Elliptical scanning	Off
Reordering	Linear
Asymmetric echo	Off
Contrasts	6
Flow comp. 1	No
Readout mode	Monopolar
Optimization	Min. TE TR
Multi-slice mode	Sequential
Echo spacing	14.4 ms
Sequence type	Gre
Bandwidth 1	300 Hz/Px
Bandwidth 2	1953 Hz/Px
Bandwidth 3	1953 Hz/Px
Bandwidth 4	1953 Hz/Px
Bandwidth 5	1953 Hz/Px
Bandwidth 6	1953 Hz/Px

Sequence - Part 2

Define	Segments
--------	----------

Sequence - Part 2

Segments	160
RF pulse type	Fast
Gradient mode	Performance*
Excitation	Non-sel.
Flip angle mode	Constant
RF spoiling	On
Phase Enc. Rewinder	On
Cine	Off

Sequence - Special

3D Rad Acq. [L] # 0	1 n.a.
3D Rad Acq. [L] # 1	2 n.a.
3D Rad Acq. [L] # 2	2 n.a.
3D Rad Rec. [L] # 0	2 n.a.
3D Rad Rec. [L] # 1	4 n.a.
3D Rad Rec. [L] # 2	2 n.a.
3D Rad Rec. [L] # 3	2 n.a.
3D Rad Rec. [L] # 4	75 n.a.
3D Rad Rec. [D] # 0	1.0 n.a.
3D Rad Rec. [D] # 1	10.0 n.a.
Prep Pulses [L] # 0	0 n.a.
LIBRE Pulse [L] # 0	1 n.a.
LIBRE Pulse [L] # 1	0 n.a.
LIBRE Pulse [L] # 2	300 n.a.
LIBRE Pulse [L] # 3	200 n.a.
LIBRE Pulse [L] # 4	0 n.a.
LIBRE Pulse [L] # 5	0 n.a.
LIBRE Pulse [L] # 6	260 n.a.
Self-Navigation [L] # 0	2 n.a.
Self-Navigation [L] # 1	2 n.a.
Self-Navigation [L] # 2	1 n.a.
Self-Navigation [L] # 3	0 n.a.
Self-Navigation [L] # 4	2 n.a.
Self-Navigation [L] # 5	2 n.a.
Self-Navigation [L] # 6	80 n.a.
Self-Navigation [L] # 7	20 n.a.
Self-Navigation [L] # 8	20 n.a.
Self-Navigation [L] # 9	1 n.a.
Self-Gating [L] 0	2 n.a.

Sequence - Assistant

Mode	Off
Allowed delay	0 s

\\QIS_LAB\BME Lab\2026\Group 4\se_multiecho_TR1500ms

TA: 3:18 PM: FIX Voxel size: 2.0×2.0×4.0 mmRel. SNR: 1.00 :

Properties

Prio recon	Off
Load images to viewer	On
Inline movie	Off
Auto store images	On
Load images to stamp segments	Off
Load images to graphic segments	Off
Auto open inline display	Off
Auto close inline display	Off
Start measurement without further preparation	Off
Wait for user to start	Off
Start measurements	Single measurement

Routine

Slice group	1
Slices	1
Dist. factor	0 %
Position	L1.7 P62.8 H24.3 mm
Orientation	Coronal
Phase enc. dir.	R >> L
AutoAlign	Head > Brain
FoV read	256 mm
FoV phase	100.0 %
Slice thickness	4 mm
TR	1500 ms
TE 1	15 ms
TE 2	30 ms
TE 3	45 ms
TE 4	60 ms
TE 5	75 ms
TE 6	90 ms
TE 7	105 ms
TE 8	120 ms
TE 9	135 ms
TE 10	150 ms
TE 11	165 ms
TE 12	180 ms
TE 13	195 ms
TE 14	210 ms
TE 15	225 ms
TE 16	240 ms
TE 17	255 ms
TE 18	270 ms
TE 19	285 ms
TE 20	300 ms
TE 21	315 ms
TE 22	330 ms
TE 23	345 ms
TE 24	360 ms
TE 25	375 ms
TE 26	390 ms
TE 27	405 ms
TE 28	420 ms
TE 29	435 ms
TE 30	450 ms
TE 31	465 ms
TE 32	480 ms
Averages	1
Filter	None
Coil elements	BO1-3;SP2-4

Contrast - Common

TR	1500 ms
TE 1	15 ms
TE 2	30 ms
TE 3	45 ms
TE 4	60 ms
TE 5	75 ms
TE 6	90 ms
TE 7	105 ms
TE 8	120 ms
TE 9	135 ms
TE 10	150 ms
TE 11	165 ms
TE 12	180 ms
TE 13	195 ms
TE 14	210 ms
TE 15	225 ms
TE 16	240 ms
TE 17	255 ms
TE 18	270 ms
TE 19	285 ms
TE 20	300 ms
TE 21	315 ms
TE 22	330 ms
TE 23	345 ms
TE 24	360 ms
TE 25	375 ms
TE 26	390 ms
TE 27	405 ms
TE 28	420 ms
TE 29	435 ms
TE 30	450 ms
TE 31	465 ms
TE 32	480 ms

Contrast - Dynamic

Averages	1
Measurements	1

Resolution - Common

FoV read	256 mm
FoV phase	100.0 %
Slice thickness	4 mm
Base resolution	128
Phase resolution	100 %

Resolution - Filter Image

Image Filter	Off
Distortion Corr.	Off
Prescan Normalize	Off
Normalize	Off
B1 filter	Off

Resolution - Filter Rawdata

Raw filter	Off
Elliptical filter	Off

Geometry - Common

Slice group	1
Slices	1
Dist. factor	0 %

Geometry - Common

Position	L1.7 P62.8 H24.3 mm
Orientation	Coronal
Phase enc. dir.	R >> L
FoV read	256 mm
FoV phase	100.0 %
Slice thickness	4 mm
TR	1500 ms
Series	Ascending

Geometry - AutoAlign

Slice group	1
Position	L1.7 P62.8 H24.3 mm
Orientation	Coronal
Phase enc. dir.	R >> L
AutoAlign	Head > Brain
Initial Position	L1.7 P62.8 F11.7
L	1.7 mm
P	62.8 mm
H	11.7 mm
Initial Rotation	0.00 deg
Initial Orientation	Coronal

Geometry - Tim Planning Suite

Set-n-Go Protocol	Off
Table position	H
Table position	36 mm
Inline Composing	Off

System - Miscellaneous

Positioning mode	FIX
Table position	H
Table position	36 mm
MSMA	S - C - T
Sagittal	R >> L
Coronal	A >> P
Transversal	F >> H
Coil Combine Mode	Adaptive Combine
Save uncombined	Off
Matrix Optimization	Off
Coil Focus	Flat
AutoAlign	Head > Brain
Coil Select Mode	Default

System - Adjustments

B0 Shim mode	Standard
B1 Shim mode	TrueForm
Adjust with body coil	Off
Confirm freq. adjustment	Off
Assume Dominant Fat	Off
Assume Silicone	Off
Adjustment Tolerance	Auto

System - Adjust Volume

! Position	L1.5 P64.8 H15.0 mm
! Orientation	S > T1.0 > C0.2
! Rotation	-0.94 deg
! A >> P	44 mm
! F >> H	137 mm
! R >> L	87 mm
Reset	Off

System - pTx Volumes

B1 Shim mode	TrueForm
--------------	----------

System - Tx/Rx

Frequency 1H	123.229536 MHz
Correction factor	1
Gain	High
Img. Scale Cor.	1.000
Reset	Off
? Ref. amplitude 1H	0.000 V

Sequence - Part 1

Contrasts	32
Bandwidth	130.208333 Hz/Px

Sequence - Assistant

Mode	Off
------	-----

\\QIS_LAB\BME Lab\2026\Group 4\se_multiecho_TR5000ms

TA: 10:45 PM: FIX Voxel size: 2.0×2.0×4.0 mmRel. SNR: 1.00 :

Properties

Prio recon	Off
Load images to viewer	On
Inline movie	Off
Auto store images	On
Load images to stamp segments	Off
Load images to graphic segments	Off
Auto open inline display	Off
Auto close inline display	Off
Start measurement without further preparation	Off
Wait for user to start	Off
Start measurements	Single measurement

Routine

Slice group	1
Slices	1
Dist. factor	0 %
Position	L1.7 P62.8 H24.3 mm
Orientation	Coronal
Phase enc. dir.	R >> L
AutoAlign	Head > Brain
FoV read	256 mm
FoV phase	100.0 %
Slice thickness	4 mm
TR	5000 ms
TE 1	15 ms
TE 2	30 ms
TE 3	45 ms
TE 4	60 ms
TE 5	75 ms
TE 6	90 ms
TE 7	105 ms
TE 8	120 ms
TE 9	135 ms
TE 10	150 ms
TE 11	165 ms
TE 12	180 ms
TE 13	195 ms
TE 14	210 ms
TE 15	225 ms
TE 16	240 ms
TE 17	255 ms
TE 18	270 ms
TE 19	285 ms
TE 20	300 ms
TE 21	315 ms
TE 22	330 ms
TE 23	345 ms
TE 24	360 ms
TE 25	375 ms
TE 26	390 ms
TE 27	405 ms
TE 28	420 ms
TE 29	435 ms
TE 30	450 ms
TE 31	465 ms
TE 32	480 ms
Averages	1
Filter	None
Coil elements	BO1-3;SP2-4

Contrast - Common

TR	5000 ms
TE 1	15 ms
TE 2	30 ms
TE 3	45 ms
TE 4	60 ms
TE 5	75 ms
TE 6	90 ms
TE 7	105 ms
TE 8	120 ms
TE 9	135 ms
TE 10	150 ms
TE 11	165 ms
TE 12	180 ms
TE 13	195 ms
TE 14	210 ms
TE 15	225 ms
TE 16	240 ms
TE 17	255 ms
TE 18	270 ms
TE 19	285 ms
TE 20	300 ms
TE 21	315 ms
TE 22	330 ms
TE 23	345 ms
TE 24	360 ms
TE 25	375 ms
TE 26	390 ms
TE 27	405 ms
TE 28	420 ms
TE 29	435 ms
TE 30	450 ms
TE 31	465 ms
TE 32	480 ms

Contrast - Dynamic

Averages	1
Measurements	1

Resolution - Common

FoV read	256 mm
FoV phase	100.0 %
Slice thickness	4 mm
Base resolution	128
Phase resolution	100 %

Resolution - Filter Image

Image Filter	Off
Distortion Corr.	Off
Prescan Normalize	Off
Normalize	Off
B1 filter	Off

Resolution - Filter Rawdata

Raw filter	Off
Elliptical filter	Off

Geometry - Common

Slice group	1
Slices	1
Dist. factor	0 %

Geometry - Common

Position	L1.7 P62.8 H24.3 mm
Orientation	Coronal
Phase enc. dir.	R >> L
FoV read	256 mm
FoV phase	100.0 %
Slice thickness	4 mm
TR	5000 ms
Series	Ascending

Geometry - AutoAlign

Slice group	1
Position	L1.7 P62.8 H24.3 mm
Orientation	Coronal
Phase enc. dir.	R >> L
AutoAlign	Head > Brain
Initial Position	L1.7 P62.8 F11.7
L	1.7 mm
P	62.8 mm
H	11.7 mm
Initial Rotation	0.00 deg
Initial Orientation	Coronal

Geometry - Tim Planning Suite

Set-n-Go Protocol	Off
Table position	H
Table position	36 mm
Inline Composing	Off

System - Miscellaneous

Positioning mode	FIX
Table position	H
Table position	36 mm
MSMA	S - C - T
Sagittal	R >> L
Coronal	A >> P
Transversal	F >> H
Coil Combine Mode	Adaptive Combine
Save uncombined	Off
Matrix Optimization	Off
Coil Focus	Flat
AutoAlign	Head > Brain
Coil Select Mode	Default

System - Adjustments

B0 Shim mode	Standard
B1 Shim mode	TrueForm
Adjust with body coil	Off
Confirm freq. adjustment	Off
Assume Dominant Fat	Off
Assume Silicone	Off
Adjustment Tolerance	Auto

System - Adjust Volume

! Position	L1.5 P64.8 H15.0 mm
! Orientation	S > T1.0 > C0.2
! Rotation	-0.94 deg
! A >> P	44 mm
! F >> H	137 mm
! R >> L	87 mm
Reset	Off

System - pTx Volumes

B1 Shim mode	TrueForm
--------------	----------

System - Tx/Rx

Frequency 1H	123.229536 MHz
Correction factor	1
Gain	High
Img. Scale Cor.	1.000
Reset	Off
? Ref. amplitude 1H	0.000 V

Sequence - Part 1

Contrasts	32
Bandwidth	130.208333 Hz/Px

Sequence - Assistant

Mode	Off
------	-----